





Recombinant Human Mucin-13(MUC13),partial (Active)

Product Code	CSB-YP887973HU
Abbreviation	Recombinant Human MUC13 protein, partial (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H3R2
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 μm filtered 20 mM Tris-HCl, 0.5 M NaCl, 6% Trehalose, pH 8.0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Biological Activity	Measured by its binding ability in a functional ELISA. Immobilized Human MUC13 at 2 μg/mL can bind Anti-MUC13 recombinant antibody (CSB-RA887973MA1HU). The EC50 is 2.592-2.923 ng/mL.
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	TNQGNSADAVTTTETATSGPTVAAADTTETNFPETASTTANTPSFPTATSPAP PIISTHSSSTIPTPAPPIISTHSSSTIPIPTAADSESTTNVNSLATSDIITASSPNDG LITMVPSETQSNEMSPPTEDNQSSGPPTGTALLETSTLNSTGPSNPCQDDPC ADNSLCVKLHNTSFCLCLEGYYNSSTCKKKGKVFPGKISVTVSETFDPEEKHS MAYQDLHSEITSLFKDVFVGTSVYGGQTVILTVSTLSRSEMRAADKDFVNVTVTI LAETTSNEKTVTEKINKAIRSSSSNFLNYDLTLRCDYYGCNQTADDCLNGLAC DCKSDLQRPNPQSPFCVASSLKCPDACNAQHKQCLIKSSGGAPECACVPGYQ EDANGNCQKCAFGYSGLDCKDKFQL
Source	yeast
Target Names	MUC13
Expression Region	19-421aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	44.6 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Activity

Measured by its binding ability in a functional ELISA. Immobilized Human MUC13 at 2 µg/ml can bind Anti-MUC13 recombinant antibody (CSB-RA887973MA1HU). The EC₅₀ is 2.592-2.923 ng/mL.

Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.