




Recombinant Mesomycoplasma hyopneumoniae 46 kDa surface antigen (p46), Biotinylated

Product Code	CSB-EP313811MWK-B
Abbreviation	Recombinant Mesomycoplasma hyopneumoniae 46 kDa surface antigen protein, Biotinylated
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C0J7
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Mesomycoplasma hyopneumoniae 46 kDa surface antigen (p46)
Immunogen Species	Mesomycoplasma hyopneumoniae (strain 232) (Mycoplasma hyopneumoniae)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	CGQTESGSTSDSKPQAETLKHKVSNDISIRIALTDPDNPRWISAQKDIISYVDET EAATSTITKNQDAQNNWLTQQANLSPAPKGFIIAPENGSGVGTAVENTIADKGIPI VAYDRLITGSDKYDWYVSFDNEKVGELQGLSLAAGLLGKEDGAFDSIDQMNE YLKSHMPQETISFYTIAGSQDDNNSQYFYNGAMKVLKELMKNSQNKIIDLSPEG ENAVYVPGWNYGTAGQRIQSFLTINKDPAGGNKIKAVGSKPASIFKGFLAPND GMAEQAITKLKLEGFDTQKIFVTGQDYNDKAKTFIKDGDQNMITYKPKVLGKV AVEVLRVLIAKKNKASRSEVENELKAKLPNISFKYDNQTYKVQGKNINTILVSPVI VTKANVDNPDA
Research Area	Others
Source	E.coli
Target Names	p46
Expression Region	28-416aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal MBP-tagged and C-terminal 6xHis-Avi-tagged
Mol. Weight	90.4 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.