




Recombinant Human Williams-Beuren syndrome chromosomal region 27 protein (WBSCR27)

Product Code	CSB-EP839810HU
Abbreviation	Recombinant Human WBSCR27 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N6F8
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Williams-Beuren syndrome chromosomal region 27 protein (WBSCR27)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MAQEEGGSLPEVRARVRAAHGIPDLAQKLHFYDRWAPDYDQDVATLLYRAPR LAVDCLTQALPGPPHSALILDVACGTGLVAAELRAPGFLQLHGVGDGSPGMLEQ AQAPGLYQRLSLCTLGQEPLPSPEGTFDAVLIVGALSDDGQVPCNAIPELHVTKP GGLVCLTTRTNSSNLQYKEALEATLDRLEQAGMWEGLVAWPVDRLWTAGSW LPPSWRWYPASLPRMASSPALSTCTESGRRPRLRK
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	WBSCR27
Expression Region	1-245aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	33.4 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.