




Recombinant Mouse Voltage-dependent L-type calcium channel subunit alpha-1D (Cacna1d), partial

Product Code	CSB-EP860754MO
Abbreviation	Recombinant Mouse Cacna1d protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99246
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Mouse Voltage-dependent L-type calcium channel subunit alpha-1D (Cacna1d), partial
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	EQGLVGKYPANKTTIALQAGLRTLHDIGPEIRRAISCDLQDDEPEDSKPEEEDV FKRNGALLGNHVNHNHNSDRRDSLQQTNTTTHRPLHVQRPSMPPASDTEKPLFP PAGNSGCHNHHNHNHNSIGKQAPTSTNANLNNANMSKAAHGKPPSIGNLEHVSE NGHYSCKHDRELQRRSSIKRTRYETIYIRSESGDEQFPTICREDPEIHGYFRD PRCLGEQEYFSSEECEDDSSPTWSRQNYNYNRYPGSSMDFERPRGYHH PQGLEDSDSPTGYDSRRSPRRRLLPPTPPSHRRSSFNFECLRRQSSQDDVL PSPALPHRAALPLHLMQQQIMAVAGLDSSKAQKYSPPSHSTRSWATPPATPPY RDWSPCYTPLIQVDRSESMDQVNGSLPSLHRSSWYTDEPDISYRTFTPASLTV PSSFRNKNSDKQRSADSLVEAVLISEGLGRYARDPKFVSATKHEIADACDLTID EMESAASTLLNGSVCPRANGDMGPISHRQ
Research Area	Cancer
Source	E.coli
Target Names	Cacna1d
Expression Region	1644-2144aa
Tag Info	C-terminal GFP-tagged
Mol. Weight	83.7 kDa
Protein Length	Partial



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.