




Recombinant *Macaca fascicularis* Parathyroid hormone/parathyroid hormone-related peptide receptor (PTH1R)), partial

Product Code	CSB-EP018988MOV3
Abbreviation	Recombinant <i>Cynomolgus</i> monkey PTH1R protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0A2K5V724
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant <i>Macaca fascicularis</i> Parathyroid hormone/parathyroid hormone-related peptide receptor (PTH1R)), partial
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	DADDVMTKEEQIFLLHRAQAQCEKRLKEVLQRPANIMESDKGWTSTSTSGKP RKDKPSGKLYPESEEDKEAPTGSRYRGRPCLPEWDHILCWPLGAPGEVVAVP CPDYIYDFNHKGHAYRRCDRNGSWELVPGHNRTWANYSECVKFLTNETRER EVFDRLG
Research Area	Others
Source	<i>E.coli</i>
Target Names	PTH1R
Expression Region	27-188aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	25.6 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.