




Recombinant Mouse Double-stranded RNA-specific adenosine deaminase (Adar), partial

Product Code	CSB-EP859990MO
Abbreviation	Recombinant Mouse Adar protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99MU3
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Mouse Double-stranded RNA-specific adenosine deaminase (Adar), partial
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	SLGTGNRCVKGDSLKGETVNDCHAEIISRGGFIRFLYSELMKYNHHTAKNSI FELARGGEKLQIKKTVSFHLYISTAPCGDGALFDKSCSDRAVESTESRHYPVFE NPKQGKLRRTKVENGEGTIPVESSDIVPTWDGIRLGERLRTMSCSDKILRWNVL GLQGALLTHFLQPVYLKSVTLGYLFSQGHLTRAIACCRVTRDGKAFEDGLRYPFI VNHPKVGRVSVYDSKRQSGKTKETSVNWC MADGYDLEILDGTRGTVDGPGK ELSRVSKKNIFLQFKKLCSEFRARRDLLQLSYGEAKKAARDYDLAKNYFKKSLRD MGYGNWISKPQEEKNF
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	Adar
Expression Region	835-1170aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	44.9 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.