




Recombinant Human Programmed cell death protein 4 (PDCD4)

Product Code	CSB-EP703610HUc7
Abbreviation	Recombinant Human PDCD4 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53EL6
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Programmed cell death protein 4 (PDCD4)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	MDVENEQILNVNPADPDNLSDSLFSGDEENAGTEEIKNEINGNWISASSINEAR INAKAKRRLRKNSSRDSGRGDSVSDSGSDALRSGTLVPTSPKGRLLDRRSRS GKGRGLPKKGGAGGKGVWGTPGQVYDVEEVDVKDPNYDDDQENCVYETVV LPLDERAFEKTLTPIIQEYFEHGDTNEVAEMLRDLNLGEMKSGVPVLAVSLALE GKASHREMTSKLLSDLCGTVMSTTDVEKSFDKLLKDLPELALDTPRAPQLVGQ FIARAVGDGILCNTYIDSYKGTVDCVQARAALDKATVLLSMSKGGKRKDSVWG SGGGQQSVNHLVKEIDMLLKEYLLSGDISEAEHCLKELEVPHFHHELVEYEAIIIM VLESTGESTFKMILDLLKSLWKSSTITVDQMKGYERYNEIPDINLDVPHSYSV LERFVEECFQAGIISKQLRDLCP SRGRKRFVSEGDGGRLKPESY
Research Area	Cancer
Source	E.coli
Target Names	PDCD4
Expression Region	1-469aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	58.7 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.