




Recombinant Mouse CTP synthase 1 (Ctps1)

Product Code	CSB-EP006176MO
Abbreviation	Recombinant Mouse Ctps1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P70698
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Mouse CTP synthase 1 (Ctps1)
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MKYILVTGGVISGIGKGVIASSVGTILKSCGLHVTSIKIDPYINIDAGTFSPYEHGE VFVLDGGEVDLDLGNRYERFLDIRLTKDNNLTTGKIYQYVINKERKGDYLGKTV QVVPHITDAIQEWVMRQALIPVDEEDGLEPQVCVIELGGTVGDIESMPFIEAFRQ FQFKVKRENFNCNIHVSLVPQPSSTGEQKTKPTQNSVRELRLGLLSPDLVVCRC SNPLDTSVKEKISMFCHEVEPEQVICVHDVSSIYRVPLLEEQQGVVDYFLRRDL PIERQSRKMLMKWKEMADRYDRLLLETCSIALVGKYTKFSDSYASVIKALEHSAL AINHKLEIKYIDSTDLEPSTLQEEPVRVYHEAWQKLCSAHGVLVPGGFGVVRGTE GKIQAIAWARKQKKPFLGVCLGMQLAVVEFVRNVLGWQDANSTEFDPKTSHP VVIDMPEHNPGQMGGTMRGKRRTLFQTKNSVMRKLYGDDTYLEERHRHRF EVNPVLKCCLEEQGLKFVGGQDVEGERMEIVELEDHPFFVGVQYHPEFLSRPIK PSPPYFGLLLASVGRPLPHYLQKGCRLSPRDTYSDRSGSSSPDSEITELKFPSIS QD
Research Area	Signal Transduction
Source	E.coli
Target Names	Ctps1
Expression Region	1-591aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	73.6 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.