




Recombinant Human Vitamin K-dependent protein C (PROC)

Product Code	CSB-EP018742HUc7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04070
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human Vitamin K-dependent protein C (PROC)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	ANSFLEELRHSSLERECIEEICDFEEAKEIFQNVDDTLAFWSKHVDGDQCLVLP LEHPCASLCCGHGTCIDGIGSFSCDCRSWEGRFQCQREVSFLNCSLDNGGCT HYCLEEVGWRRCSAPGYKLGDDLQCHPAVKFPCGRPWKMEKKRSHLKR DTEDQEDQVDPRLIDGKMTRRGDSPWQVLLDSKKKLACGAVLIHPSWVLT AHCMDSESKLLVRLGEYDLRRWEKWELDLDIKEVFPVHPNYSKSTTDNDIALH LAQPATLSQTIVPICLPDSGLAERELNQAQETLVTGWGYHSSREKEAKRNRT FVLNFIKIPVPHNECSEVMSNMVSENMLCAGILGDRQDACEGDSGGPMVAS FHGTWFLVGLVSWGEGCGLLHNYGVYTKVSRYLWDWIHGHIRDKEAPQKSWA P
Research Area	Others
Source	E.coli
Target Names	PROC
Expression Region	43-461aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	48.3 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.