




Recombinant Pig Protocadherin beta 7 (PCDHB7), partial

Product Code	CSB-EP5745PI
Abbreviation	Recombinant Pig PCDHB7 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0A287BG69
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Pig Protocadherin beta 7 (PCDHB7), partial
Immunogen Species	Sus scrofa (Pig)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	GAEPLRYFVEEETERGTFLANLANDLGLGVGELSARAPRIGSDQNIRFLLLNP TGDLLLNEKLDREELCGPAEPCVLPFQLLLEEPLEIFRAELWVRDINDNSPVFL HREILLKILESTTAGAAFLLESAQDSAGTNSLRNYTIGPNDYFHITVHDGEEGIL YPELVMDRALDREQVPELSLTLTALDGGSPPRSGTALVHIQVLDVNDNAPEFV QSLYKVQVPENSPVGSVVTVSARDLDAGSNGEIVYAFLYATERILETFQINPA SGSLYLKAEKNYEAVQTYTLIIQAKDGGGLSGKCTVVVEVTDINDNPPELVMSS LTSSIKENSPETVVAVFRIRDRDSGNNAKMVCISIAHLFPVLKPSVENYYTLVT ERPLDRESQAEYNITISVSDLGTPRRHTQHNLTVRVADVNDNAPAFSQPAYTL RVRENNSPALHIGTVRATDADAGANAQLTYSLLPPSDPHLPLASLVAINPDTGQ LFALRALDYEALRAFDQVRAADRGAPALSSQARVRVLVADANDHAPVVLYPP PNASAPCTELLPRAAEPGYLLAKVVAVDGDAGQNAWLSYRLLQATEPGLFGV GAHSGEVRTARPLSERDAPRQRLLLLVRDHGEPPLSASVPLHVLLVDGFSQPY LPPPEAPGRAPGAGAGADEP
Research Area	Neuroscience
Source	E.coli
Target Names	PCDHB7
Expression Region	21-684aa
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	78.8 kDa
Protein Length	Partial



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.