




AMP2-2 Antibody

Product Code	CSB-PA871127XA01MAB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9SPL4
Immunogen	Recombinant Macadamia integrifolia (Macadamia nut) AMP2-2 protein (271-410aa)
Raised In	Rabbit
Species Reactivity	Macadamia integrifolia (Macadamia nut)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Macadamia integrifolia (Macadamia nut)
Target Names	AMP2-2 Antibody

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes:AMP2-2 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 18 kDa

Observed band size: 18 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.