




forkhead box protein O Antibody

Product Code	CSB-PA5622XA01
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	XP_047510659.1
Immunogen	Recombinant Pieris napi forkhead box protein O protein (281-470aa)
Raised In	Rabbit
Species Reactivity	Pieris napi
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pieris napi
Target Names	forkhead box protein O Antibody

Image



Western Blot

Positive WB detected in: recombinant protein
All lanes: forkhead box protein O Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution
Predicted band size: 35.6 kDa
Observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.