



Phospho-PSEN1 (T367) Antibody

Product Code	CSB-PA5730XA01phPI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F1S3K9
Immunogen	Peptide sequence around phosphorylation site of serine T367 from Human PSEN1.
Raised In	Rabbit
Species Reactivity	Sus scrofa (Pig)
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sus scrofa (Pig)
Target Names	Phospho-PSEN1 (T367) Antibody
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.