





DPEP3 Recombinant Monoclonal Antibody

Product Code	CSB-RA007125MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Human DPEP3 protein
Species Reactivity	Human, Macaca fascicularis
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	mIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DPEP3
Clone No.	5H3

Image



The Binding Activity of Human DPEP3 with Anti-DPEP3 recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human DPEP3 (CSB-MP007125HU) at 2 µg/mL can bind Anti-DPEP3 recombinant antibody. The EC₅₀ is 6.841-7.498 ng/mL.



The Binding Activity of Macaca fascicularis DPEP3 with Anti-DPEP3 recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Macaca fascicularis DPEP3 (CSB-MP007125MOV) at 2 µg/mL can bind Anti-DPEP3 recombinant antibody. The EC₅₀ is 7.817-8.936 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.