




CDH1 Recombinant Monoclonal Antibody

Product Code	CSB-RA005034MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Human CDH1 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer;Developmental biology;Signal transduction
Target Names	CDH1
Clone No.	10H4

Image



The Binding Activity of Human CDH1 with Anti-CDH1 recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CDH1(CSB-MP005034HU1) at 2 µg/mL can bind Anti-CDH1 recombinant antibody. The EC₅₀ is 5.151-6.358 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.