





IL1RAP Recombinant Monoclonal Antibody

Product Code	CSB-RA878844MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPH3
Immunogen	Recombinant Human IL1RAP protein
Species Reactivity	Human, Macaca fascicularis
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular;Immunology
Target Names	IL1RAP
Clone No.	2G7

Image



The Binding Activity of Macaca fascicularis IL1RAP with Anti-IL1RAP recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Macaca fascicularis IL1RAP (CSB-MP5268MOV) at 2 µg/mL can bind Anti-IL1RAP recombinant antibody. The EC₅₀ is 0.9099-1.181 ng/mL.



The Binding Activity of Human IL1RAP with Anti-IL1RAP recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human IL1RAP (CSB-MP878844HU) at 2 µg/mL can bind Anti-IL1RAP recombinant antibody. The EC₅₀ is 0.9667-1.477 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.