




IL12RB1 Recombinant Monoclonal Antibody

Product Code	CSB-RA011588MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Human IL12RB1 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL12RB1
Clone No.	7C11

Image



The Binding Activity of Human IL12RB1 with Anti-IL12RB1 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized human IL12RB1(CSB-MP011588HU3) at 2 µg/mL can bind Anti-IL12RB1 recombinant antibody. The EC₅₀ is 1.578-2.143 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.