





IL2RG Recombinant Monoclonal Antibody

Product Code	CSB-RA011651MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Human IL2RG protein
Species Reactivity	Human, Cynomolgus
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG4
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer;Immunology
Target Names	IL2RG
Clone No.	6F3

Image



The Binding Activity of Rhesus macaque IL2RG with Anti-IL2RG recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Rhesus macaque IL2RG (CSB-MP2765MOW) at 2 µg/mL can bind Anti-IL2RG recombinant antibody. The EC₅₀ is 33.73-41.04 ng/mL.



The Binding Activity of Human DLL3 with Anti-DLL3 recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human DLL3 (CSB-MP882142HU2d7) at 2 µg/mL can bind Anti-DLL3 recombinant antibody. The EC₅₀ is 3.990-4.723 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.