




IFNA2 Recombinant Monoclonal Antibody

Product Code	CSB-RA365593MA1MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01573
Immunogen	Recombinant Mouse IFNA2 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	IFNA2
Clone No.	8E8

Image



The Binding Activity of Human IFNA2with Anti-IFNA2 Recombinant Antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human IFNA2(CSB-MP360706HU) at 2 µg/mL can bind Anti-IFNA2 recombinant antibody. The EC₅₀ is 2.366-2.818 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.