




LGALS3 Recombinant Monoclonal Antibody

Product Code	CSB-RA012887MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17931
Immunogen	Recombinant Human LGALS3 protein
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience;Cell biology
Target Names	LGALS3
Clone No.	5E5

Image



Western Blot

Positive WB detected in: PC-3 whole cell lysate, MCF7 whole cell lysate, HeLa whole cell lysate, A431 whole cell lysate

All lanes: LGALS3 antibody at 1:1000

Secondary

Goat polyclonal to human IgG at 1/50000 dilution

Predicted band size: 26 kDa

Observed band size: 26 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.