




Recombinant *Xenopus laevis* Progesterone and adiponectin receptor family member VIII L homeolog (paqr8.L)-VLPs

Product Code	CSB-MP4097XBE
Abbreviation	Recombinant <i>Xenopus laevis</i> paqr8.L protein-VLPs
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q801G4
Form	Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant <i>Xenopus laevis</i> Progesterone and adiponectin receptor family member VIII L homeolog (paqr8.L)-VLPs
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Sensitivity	Not Test
Purity	VLPs
Sequence	MTTAILECISTLSISFQQLRRLPRFLEGGTTKMPLTVTDSVPRLFREPIYQTGY RPTDQDWKYYFLSLFKKHNESVNVWTHLLVALAVVLRVVAFFVEAGSLSLNVVS FPLYLYVLSSLTYLTCSILAHLLQSKSELAHYTFYFIDYVGVSTYQYGCALAHYY YTSNEAWYDKACYFFLPGAFLGWLSCVGCCYAKYCYKRPYPVMRKILQVVP AGLAYILDISPVVHRIVTCHMEDYTDKAVWLHSLQMIFFIIGAYFFSCPVPKEYF PGSCDFIGHGHQIFHVFLGLCTLSQLEALFIDYQTRQEVFSARYSSNYTLMCCA SFFLLILCSTFTAVYARRRIKEKLARKEL
Research Area	Others
Source	Mammalian cell
Target Names	paqr8.L
Expression Region	1-353aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	41.9 kDa
Protein Length	Full Length
Image	



CSB-MP4097XBE is detected by Mouse anti-6*His monoclonal antibody.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.