




Recombinant Human herpesvirus 1 Triplex capsid protein 2 (TRX2)

Product Code	CSB-MP320726HWY
Abbreviation	Recombinant Human herpesvirus 1 TRX2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P10202
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Human herpesvirus 1 Triplex capsid protein 2 (TRX2)
Immunogen Species	Human herpesvirus 1 (strain 17) (HHV-1) (Human herpes simplex virus 1)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MLADGFETDIAIPSGISRPDAAALQRCEGRVVFLPTIRRQLTLADVAHESFVSG GVSPDTLGLLLAYRRRFPVITRVLPTRIVACPLDVGLTHAGTVNLRNTSPVDL CNGDPISLVPPVFEGQATDVRLDSLDTLRFPVPLPSPLAREIVARLVARGIRDL NPSPRNPGGLPDLNVLYYNGSRLSLLADVQQLGPVNAELRSLVLNMVYSITEG TTIILTLIPRLFALSAQDGYVNALLQMQSVTREAAQLIHPEAPALMQDGERRLPL YEALVAWLTHAGQLGDTLALAPVVRVCTFDGAAVVRSGDMAPVIRYP
Research Area	Others
Source	Mammalian cell
Target Names	TRX2
Expression Region	1-318aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal hFc1-tagged
Mol. Weight	63.1 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.