




Recombinant Dog Albumin (ALB), Biotinylated

Product Code	CSB-EP001561DOg2-B
Abbreviation	Recombinant Dog ALB protein, Biotinylated
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P49822
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose.
Product Type	Recombinant Dog Serum albumin (ALB), Biotinylated
Immunogen Species	Canis familiaris(Dog)(Canis lupus familiaris)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	EAYKSEIAHRYNDLGEEHFRGLVLVAFSQYLQQCPFEDHVKLAKEVTEFAKAC AAESGANCDKSLHTLFGDKLCTVASLRDKYGDMADCCEKQEPDRNECFLAH KDDNPGFPPLVAPEPDALCAAFQDNEQLFLGKLYEIARRHPYFYAPELLYYA QQYKGVFAECCQAADKAAACLGPKIEALREKVLLSSAKERFKCASLQKFGDRAF KAWSVARLSQRFPKADFAEISKVVTDLTQVHKECCHGDLLECADDRADLAKY MCENQDSISTKLKECCDKPVLEKSQCLAEVERDELPGDLPSLAADFVEDKEVC KNYQEAKDVFLGTFLYEYARRHPEYSVSLLLRLAKEYEATLEKCCATDDPPTC YAKVLDEFKPLVDEPQNLVKTNCELFEKLGEYGFQNALLVRYTKKAPQVSTPT LVEVSRKLGKVGTKCCKKPESERMSCAEDFLSVVLNRLCVLHEKTPVSEKVTK CCSESLVNRRPCFSGLEVDETYVPKEFNAETFTFHADLCTLPEAEKQVKKQTA LVELLKHKPKATDEQLKTMGDFGAFVEKCCAAENKEGCFSEEGPKLVAAAQ AALV
Research Area	Cardiovascular
Source	E.coli
Target Names	ALB
Expression Region	25-608aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal6xHis-Avi-tagged
Mol. Weight	68.7
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.