




Recombinant Mouse Nuclear factor interleukin-3-regulated protein (Nfil3)

Product Code	CSB-YP015757MO
Abbreviation	Recombinant Mouse Nfil3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O08750
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Mouse Nuclear factor interleukin-3-regulated protein (Nfil3)
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	MQLRKMQTIKKEPAPLDPTSSSDKMLLLNSALAEVAEDLASGEDLLLNEGSMG KNKSSACRRKREFIPDEKKDAMYWEKRRKNNEAAKRSREKRRLNDLVLENKLI ALGEENATLKAELLSLKLKFLISSTAYAQEIQKLSNSTAVYFQDYQTSKAAVSS FVDEHEPAMVAGSCISVIKHSPQSSLSDVSEVSSVEHTQESPAQGGCRSPEN KFPVIKQEPVELESFAREAREERGTYSTSIYQSYMGSSFSTYSHSPPLLQVHG STSNSPRTSEADEGVVGKSSDGEDEQQVPGPIHSPVELQRVHATVVKVPEV NPSALPHKLRIKAKAMQVKVEALDSEFEGMQKLSSPADAIKRHFDLEKHGTS GMAHSSLPPFSVQVTNIQDWSLKSEHWHHKELSSKTQSSFKTGVVEVKDGGY KVSEAENLYLKQGIANLSAEVVSLKRFIATQPISASDSR
Research Area	Epigenetics and Nuclear Signaling
Source	Yeast
Target Names	Nfil3
Expression Region	1-462aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	52.7 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.