



Recombinant Pyrococcus furiosus Cell division control protein 6 homolog (cdc6)

Product Code	CSB-MP005022FHW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P81413
Product Type	Recombinant Protein
Immunogen Species	Pyrococcus furiosus (strain ATCC 43587 / DSM 3638 / JCM 8422 / Vc1)
Purity	≥85% (SDS-PAGE)
Sequence	MNEGEHQIKL DELFEKLLRA RKIFKNKDV L RHSYTPKDLP HRHEQIETLA QILVPVLRGE TPSNIFVYGK TGTGKT VTKV FVTEELKRIS EKYNIPVDVI YINCEIVDTH YRVLANIVNY FKDETGIEVP MVGWPTDEVY AKLKQVIDMK ERFVIIVLDE IDKLVKKS GD EVLYSLTRIN TELKRAKVS V IGISNDLKFK EYLDPRVLSS LSEEEVFP YDANQLRDIL TQRAEEAFYP GVLDEGVIPL CAALAAREHG DARKALDLLR VAGEIAEREG ASKVTEKHVW KAQEKIEQDM MEEVIKTLPL QSKVLLYAI V LLDENGDLPA NTGDVYAVYR ELCEYIDLEP LTQRRISDLI NELDMLGIIN AKVVS KGRYG RTKEIRLNVT SYKIRNV LRY DYSIQPL LTI SLKSEQRRLI
Source	Mammalian cell
Target Names	cdc6
Protein Names	Recommended name: Cell division control protein 6 homolog Short name= CDC6 homolog
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.