



Recombinant Mouse Nucleolar RNA helicase 2 (Ddx21)

Product Code	CSB-EP881308MO(A4)c7
Abbreviation	Recombinant Mouse Ddx21 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9JIK5
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Mouse Nucleolar RNA helicase 2 (Ddx21)
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MPGKLRSKAKLGSAGAEESMETLPKPSEKKTRKEKTKSKTEEATEGMEEAVS SKAKKTNKKGPSEDDVDPPKSRKAKKQEEEPQDDTASTSKTSKKKKEPLEKQ ADSETKEIITEEPSEEEADMPKPKMKKGKEANGDAGEKSPKLKNGLSQPSEE EADIPKPKMKKGKEANGDAGEKSPKLKNGLSQPSEEEVDIPKPKMKKGKE ASGDAGEKSPRLKDGLSQPSEPKSNSSDAPGEESSETEKEIPVEQKEGAFS NFPISEETVKKLKARGVNFLFPIQAKTFHHVYSGKDLIAQARTGTGKTFSAIPLI EKLQGGGLQERKRGRAPQVLVLAPTRELANQVSKDFSDITKKLSVACFYGGTPY GGQIERMRSGIDILVGTPGRIKDHLQNGKLDLTKLKHVVLDEVDQMLDMGFAD QVEEILCVAYKKDSEDNPQTLLFSATCPHWVFNVAKKYMKSTYEQVDLIGKKT QKAAITVEHLAIKCHWTERAAVIGDVIRVYSGHQGRTIIFCETKKDAQELSQNTC IKQDAQSLHGDIPQKQREITLKGFRNGNFGVLVATNVAARGLDIPEVDLVVQSC PPKDVESYIHRSGRTGRAGRTGVCICFYQNKEEYQLAQVEQKAGIKFKRIGVP SATEIHKASSKDAIRLLDSVPPTAISHFQKSAEKLIEEKGAVEALAAALAHISGATS VDQRSLINSQAGFVTMILRCSIEMPNISYAWKELKEQLGESIDAKVKGMVFLKG KLGVCFDV RTEAVTEIQEKWHDSSRRWQLTVATEQPELEGPPDGYRGRMGQR DGSRGAFRGQRGGSRNFRGQGQRGGSRNFRGQRPGGGNRRGQKRSFSKAF GQ
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	Ddx21
Expression Region	1-851aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at



4°C for up to one week.

Tag Info


C-terminal 6xHis-tagged

Mol. Weight

100.5 kDa

Protein Length

Full Length

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.