



Recombinant Monkeypox virus Soluble TNF receptor II (OPG002)

Product Code	CSB-EP286411MHV-B
Uniprot No.	P0DTN0
Immunogen Species	Monkeypox virus
Sequence	MRSVLYSYILFLSIIINGRDIAPHAPSNGKCKDNEYRSRNLCCCLSCPPGTYASR LCDSKTNTQCTPCGSDTFTSHNNHLQAACLSNRCDSNQVETRSCNTTHNRI CECSPGYCYLLKGASGCRTCISKTKCGIGYGVSGYTSTGDVICSPCGPGTYSH TVSSTDKCEPVVTSNTFNIDVEINLYPVNDTSTRTTTTGLSEISITSELTITMN HKDCDPVFRAEYFVSLNNVATSGFFTGENRYQNTSKICTLNFEIKCNNKDSSS KQLTKTKNDTIMPHSETVTLVGDCLSSVDIYILYSNTNTQDYENDTISYHMGNV LDVNSHMPASCDIHKLITNSQNPTHL
Source	E.coli
Target Names	OPG002
Expression Region	20-349aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.