



Recombinant Monkeypox virus Host range factor p28 (OPG021)

Product Code	CSB-EP027267MHV
Uniprot No.	Q8V571
Immunogen Species	Monkeypox virus (strain Zaire-96-I-16) (MPX)
Sequence	MEFDPVKINTSSIDHVTLQYIDEPNDIRLTVCIIRNINNITYYINITKINPHLANRFR AWKKRIAGRDMYMTNLSRDTGIQQSKLTETIRNCQKNKNIYGLYIHYNLVINVID WITDVIVQSILRGLVNWYIANNTYTPNTPNNTTTISELDIIKILDKYEDMYRVSKE KEGCICYEVVYSKRENDRYFGLLDSCNHIFCITCINIWHRRTRRETGALDNCPIC RTRFRKITMSKFYKLVN
Source	E.coli
Target Names	OPG021
Expression Region	1-242aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.