



Recombinant Phleum pratense Profilin-1 (PRO1)

Product Code	CSB-YP017824EUQ
Abbreviation	Recombinant Phleum pratense PRO1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35079
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Immunogen Species	Phleum pratense (Common timothy)
Sensitivity	Not Test
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	SWQTYVDEHLMCEIEGHHLASAAILGHDGTVWAQSADFPQFKPEEITGIMKDF DEPGHLAPTGMFVAGAKYMWIQGEPRVIRGKKGAGGITIKKTGQALVVGIYD EPMTPGQCNCMVVERLGDYLVLEQGM
Research Area	Others
Source	Yeast
Target Names	PRO1
Expression Region	2-131aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	15.6 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°C. Solubilize for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.