



Recombinant Bovine Seizure 6-like protein 2 (SEZ6L2), partial


Product Code	CSB-YP643717BO
Abbreviation	Recombinant Bovine SEZ6L2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q29RN8
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Immunogen Species	Bos taurus (Bovine)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	LPLKEEEALPEPGSETPTVASEALAE LLHGALLRKGP EMGYLPGSDPDPTLAT PPVGQTIAAPFLPRATEPGTGPLTTAVTPKGGRGAGPTAPELLTPPPGTTAPPL PGPASPGPPLGPEGGEEETTTTITTTT VTTT VTSPLC NNNISEGEGHVESPD LGSSTSRTLGLLDCTYSIHVYPGYGIEIQVQMLNLSREEELLVLAGGGSPGLAP RLLANSSMLGEGQVLRSPNRLLLHFQSPRVPRGAGFRIHYQAYLLSCGFPPQ PAHGDVSVTDLHPGGTATFHCD SGYQLQGEETLVCLNGTRPAWSSEPPSCM ASCGGTIHNATLGRIVSPEPGGAAGPNLTCRWVIEAAEGRRLHLHFERSLDE DNDRLMVRSGGSPLSPVIYDSDMDDVPERGLISDAQSLYVELLSETPANPLLL SLRFEAFEEDRCFAPFLAHGNVTTTDPEYRPGALATFSCLPGYALEPPGPPNVI ECVDPTEPHWNDTEPACKAMCGGELSEPAGVVLSPDWPQSYSPGQDCVWG LHVQEEKRILLQVEILNVREGDMLTLFDGDGPSARVLAQLRGPQPRRLLSSG SDLT LQFQAPP GPPNPGLGQG FVLHFKEVPRNDTCPELPPPEWGWRTASHG DLIRGTVLTYQCEPGYELLGSDILTCQWDL SWSAAPPACQKIMTCADPGEITN GHRTTSDAGFPVGSHVQYRCLPGYSLEGA AVLTCYSRDTGTPKWSDRVPKC ALKYEPCLNPGVPENGYQTLYKH HYQAGESLRFFCYEGFELIGEVTITCVPGH PSQWTSQPPLCKVTQTDDPSRQLEGGN
Research Area	Neuroscience
Source	Yeast
Target Names	SEZ6L2
Expression Region	28-844aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged



Mol. Weight 89.1 kDa

Protein Length Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.