



Recombinant Rat Sphingosine 1-phosphate receptor 2 (S1pr2)-VLPs


Product Code	CSB-MP020651RA
Abbreviation	Recombinant Rat S1pr2 protein-VLPs
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47752
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Rat Sphingosine 1-phosphate receptor 2 (S1pr2)-VLPs
Immunogen Species	Rattus norvegicus (Rat)
Sensitivity	Not Test
Purity	The purity information is not available for VLPs proteins.
Sequence	MGGLYSEYLNPEKVQEHYNYTKETLDMQETPSRKVASAFIIILCCAIVVENLLVL IAVARN SKFHSAMYFLGNLAASDLLAGVAFVANTLLSGPVTLSLTPLQWFARE GSAFITLSASVFSLLAIAIERQVAIAKVKLYGSDKSCRMLMLIGASWLISLILGGLP ILGWNCLDHLEACSTVLPLYAKHYVLCVVTIFSVILLAIVALYVRIYFVVRSSHAD VAGPQTLALLKTVTIVLGVFIICWLPAFSILLDDSTCPVRACPVLYKAHYFFAFAT LNSLLNPVIYTWRSRDLRREVLRLPLLCWRQKGKATGRRGGNPGHRLLPLRSS SSLERGLHMPTSPTFLEGNTVV
Research Area	Epigenetics and Nuclear Signaling
Source	Mammalian cell
Target Names	S1pr2
Expression Region	1-352aa
Notes	The VLPs are expressed from human 293 cells (HEK293). Mix the sample gently by repeatedly pipetting it up and down. Do not vortex. Repeated freezing and thawing is not recommended. Store the protein at -20°C/-80°C upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity. The immunization strategy should be optimized (antigen dose, regimen and adjuvant).
Tag Info	C-terminal 10xHis-tagged(This tag can be tested only under denaturing conditions)
Mol. Weight	40.4 kDa



Protein Length

Full Length

Image



CSB-MP020651RA is detected by Mouse anti-6*His monoclonal antibody.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°C. Solubilize for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.