



Recombinant Rat Transient receptor potential cation channel subfamily V member 2 (Trpv2)-VLPs

| | |
|--------------------------|---|
| Product Code | CSB-MP893416RA |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9WUD2 |
| Form | Lyophilized powder |
| Storage Buffer | Lyophilized from PBS, 6% Trehalose, pH 7.4. |
| Product Type | Recombinant Rat Transient receptor potential cation channel subfamily V member 2 (Trpv2)-VLPs |
| Immunogen Species | Rattus norvegicus (Rat) |
| Sequence | <p>MTSASSPPAFRLETSDGDEEGNAEVNKGKQEPPPMESPFQREDRNSSPQIKV NLNFIKRPPKNTSAPSQQEPDRFDRDRLFSVVSARGVPEELTGLLEYLRWNSKY LTDSAYTEGSTGKTCLMKAVLNLQDGVNACIMPLLQIDKDSGNPKLLVNAQCT DEFYQGHSAHIAIEKRSLQCVKLLVENGADVHLRACGRFFQKHQGTCFYFGE LPLSLAACTKQWDVVTYLLENPHQPASLEATDSLGNLTVLHALVMIADNSPENS ALVIHMYDGLLQMGARLCPTVQLEEISNHQGLTPLKLAKEGKIEIFRHILQREF SGPYQPLSRKFTWCYGPVRSVLYDLSSVDSWEKNSVLEIIAFHCKSPNRHR MVLVLEPLNKLKQEKWDRLVSRFFNFACYLVYMFIFTVVAYHQPSLDQPAIPSS KATFGESMLLLGHILLLGGIYLLLGQLWYFWRRLFIWISFMDSYFEILFLLQAL LTVLSQVLRFMETEWYPLLVLVSLVGLWLNLLYYTRGFQHTGIYSVMIQKVILR DLLRFLVYLVFLFGFAVALVSLSREARSPKAPEDNNSTVTEQPTVGQEEEPAP YRSILDASLELFKFTIGMGELAFQEQLRFRGVLLLLLAYVLLTYVLLNMLIALM SETVNHVADNSWSIWKLQKAISVLEMENGYWWCRRKKHREGRLKVGTRGD GTPDERWCFRVEEVNWAWEKTLPTLSEDPSPGPGITGNKKNPTSKPGKNSA SEEDHLPLQVLQSP</p> |
| Research Area | Neuroscience |
| Source | Mammalian cell |
| Target Names | Trpv2 |
| Expression Region | 1-761aa |
| Notes | The VLPs are expressed from human 293 cells (HEK293). Mix the sample gently by repeatedly pipetting it up and down. Do not vortex. Repeated freezing and thawing is not recommended. Store the protein at -20°C/-80°C upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity. The immunization strategy should be optimized (antigen dose, regimen and adjuvant). |

**Tag Info**

C-terminal 10xHis-tagged

If you have specified tag type, please tell us and we will check if it's possible to develop.

Protein Length

Full Length

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°C. Solubilize for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.