




# DM1&DM4 Monoclonal Antibody

<b>Product Code</b>	CSB-MA989811I2m
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	DM4-BSA
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Homo sapiens (Human)
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Monoclonal Antibody
<b>Research Area</b>	Others
<b>Target Names</b>	DM1&DM4
<b>Clone No.</b>	12H10E10

## Image



The Binding Activity of Trastuzumab emtansine (Ado-Trastuzumab emtansine) and Cantuzumab ravtansine (Synonyms:IMGN242) with Anti-DM1&DM4 antibody

Activity: Measured by its binding ability in a functional ELISA. Immobilized Trastuzumab emtansine (Ado-Trastuzumab emtansine) and Cantuzumab ravtansine (Synonyms:IMGN242) at 2 µg/mL can bind Anti-DM1&DM4 antibody, the EC<sub>50</sub> is 4.088 to 5.365 ng/mL and 3.215 to 4.720 ng/mL, respectively.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.