




MMAE Monoclonal Antibody

Product Code	CSB-MA096584I0m
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	MMAE-BSA
Raised In	Mouse
Species Reactivity	Homo sapiens (Human)
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG1
Clonality	Monoclonal
Product Type	Monoclonal Antibody
Research Area	Others
Target Names	MMAE
Clone No.	40F3D10

Image



The Binding Activity of Disitamab vedotin (Synonyms: RC48) and Depatuxizumab mafodotin (Synonyms: ABT-414) with Anti-MMAE antibody

Activity: Measured by its binding ability in a functional ELISA. Immobilized Disitamab vedotin (Synonyms: RC48) and Depatuxizumab mafodotin (Synonyms: ABT-414) at 2 µg/mL can bind Anti-MMAE antibody, the EC₅₀ is 25.14 to 32.78 and no.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.