




Recombinant Rat Glutamate receptor ionotropic, NMDA 1 (Grin1), partial-VLPs, Fluorescent

Product Code	CSB-MP009911RA2f4
Abbreviation	Recombinant Rat Grin1 protein, partial, Fluorescent-VLPs
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P35439
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Rat Glutamate receptor ionotropic, NMDA 1 (Grin1), partial-VLPs, Fluorescent
Immunogen Species	Rattus norvegicus (Rat)
Sensitivity	Not Test
Purity	The purity information is not available for VLPs proteins.
Sequence	TRLKIVTIHQEPFVYVKPTMSDGTCKEEFTVNGDPVKKVICTGPNDTSPGSPRH TVPQCCYGFCIDLLIKLARTMNFTYEVHLVADGKFGTQERVNNSNKKKEWNGM MGELLSGQADMIVAPLTINNERAQYIEFSKPFKYQGLTILVKKEIPRSTLDSFMQ PFQSTLWLLVGLSVHVAVMLYLLDRFSPFGRFKVNSEEEEEEDALTSSAMWF SWGVLLNSGIGEGAPRSFSARILGMVWAGFAMIIVASYTANLAAFLVLD RPEER ITGINDPRLRNPSDKFIYATVKQSSVDIYFRRQVELSTMYRHMEKHNYESAAEA IQAVRDNKLHAFIWDSAVLEFEASQKCDLVTTGELFFRSGFGIGMRKDSPWKQ NVLSILKSHENGFMEDLDKTWVRYQECDS
Research Area	Others
Source	Mammalian cell
Target Names	Grin1
Expression Region	396-800aa
Notes	The VLPs are expressed from human 293 cells (HEK293). Mix the sample gently by repeatedly pipetting it up and down. Do not vortex. Repeated freezing and thawing is not recommended. Store the protein at -20°C/-80°C upon receiving it, and ensure to avoid repeated freezing and thawing, otherwise, it will affect the protein activity. The immunization strategy should be optimized (antigen dose, regimen and adjuvant).
Tag Info	C-terminal EGFP-tagged
Mol. Weight	73.3 kDa


Protein Length

Partial

Image


CSB-MP009911RA2f4 is detected by Mouse anti-GFP monoclonal antibody.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°C. Solubilize for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.