



PPP2R5C Antibody

Product Code	CSB-PA619761GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13362
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A 56 kDa regulatory subunit gamma isoform protein (1-380AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB, IHC
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	PPP2R5C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.