




SIGLEC9 Recombinant Monoclonal Antibody

Product Code	CSB-RA021303MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y336
Immunogen	Recombinant Human SIGLEC9 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular;Tags & Cell Markers?Immunology
Target Names	SIGLEC9
Clone No.	5H2

Image



The Binding Activity of Human SIGLEC9 with Anti-SIGLEC9 Recombinant Antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human SIGLEC9 (CSB-MP021303HU) at 2 µg/mL can bind Anti-SIGLEC9 recombinant antibody. The EC₅₀ is 1.033-1.429 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.