





CEACAM1 Recombinant Monoclonal Antibody

Product Code	CSB-RA005157MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13688
Immunogen	Recombinant Human CEACAM protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer;Tags & Cell Markers?Immunology;Microbiology
Target Names	CEACAM1
Clone No.	7F8

Image



The Binding Activity of Human CEACAM1 with Anti-CEACAM1 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CEACAM1 at 2 µg/mL can bind Anti-CEACAM1 recombinant antibody(CSB-RA005157MA1HU). The EC₅₀ is 3.432 - 6.079 ng/mL.



The Binding Activity of Human CEACAM1 with Anti-CEACAM1 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CEACAM1(CSB-YP005157HU) at 2 µg/mL can bind Anti-CEACAM1 recombinant antibody. The EC₅₀ is 0.3677-0.4528 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.