



Recombinant Sparus aurata Fibroblast growth factor (fgf2) (Active)

Product Code	CSB-EP6199FJX-WD
Abbreviation	Recombinant Sparus aurata FGF-2 protein (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	XP_030252529.1/A0A671VAT8
Form	Lyophilized
Storage Buffer	Lyophilized from a 0.2 μm filtered solution containing 10 mM PB, 250 mMNaCl, 5% mannitol and 0.01% Tween 80, pH7.4
Product Type	Recombinant Protein
Immunogen Species	Sparus aurata (Gilthead sea bream)
Biological Activity	Measured in a cell proliferation assay using NIH3T3 mouse embryonic fibroblast cells. The ED50 for this effect is ≤2.0 ng/mL.
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	MATGEITTLPTPEDGGSGGFPPGFSFKEPKRLYCKNGGFFLRIKSDGGVDGIR EKNDPHIKLQLQATSVGEVVIKGVCANRYLAMNRDGRFLFGARRATDECHFLER LESNNYNTYRSRKYPNMYVALKRTGQYKSGSKTGPGQKAILFLPMSAKC
Source	E.coli
Target Names	bFGF/FGF-2
Expression Region	1-155aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag free
Mol. Weight	17.1 kDa
Protein Length	Full Length
Endotoxin	≤200 EU/mg by the LAL method
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.