




Recombinant Human IL12P40&IL12P35 Heterodimer Protein (Active)

Product Code	CSB-MP6495HU-WD
Abbreviation	Recombinant Human IL12P40&IL12P35 Heterodimer protein (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29460/P29459
Form	Lyophilized
Storage Buffer	Lyophilized from a 0.2 μm filtered solution containing PBS, 5% mannitol and 0.01% Tween 80, pH 7.4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Biological Activity	Measured by its ability to induce Interferon gamma secretion by human natural killer lymphoma NK-92 cells. The ED50 for this effect is 20.5 ng/mL.
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	IWELKKDVYVVELDWYPDAPGEMVVLTCDTPEEDGITWTLTLDQSSEVLGSGKTL TIQVKEFGDAGQYTCHKGGEVLSHLLLLHKKEDGIWSTDILKDQKEPKNKFTL RCEAKNYSGRFTCWWLTTISTDLTFSVKSSRGSSDPQGVTCGAATLSAERVR GDNKEYEYSVECEQEDSACPAAEESLPIEVMVDAVHKLKYENYSSFFIRDIKP DPPKNLQLKPLKNSRQVEVSWEYPDTWSTPHSYFSLTFCVQVQGGKSKREKGD RVFTDKTSATVICRKNASISVRAQDRYSSSWSEWASVPCS&RNLPVATPDP GMFPCLHHSQNLLRAVSNMLQKARQTLEFYPTSEEIDHEDITKDKTSTVEAC LPLELTKNESCLNSRETSFITNGSCLASRKTSFMMALCLSSIYEDLKMYQVEFK TMNAKLLMDPKRQIFLDQNMLAVIDELMQALNFNSETVPQKSSLEEDPDFYKTKI KLCILLHAFRIRAVTIDRVMSYLNAS
Source	Mammalian cell
Target Names	IL-12
Expression Region	23-328aa(IL12P40)&23-219aa(IL12P35)
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag free
Mol. Weight	59.0 kDa
Protein Length	Heterodimer
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Endotoxin

≤10 EU/mg by the LAL method

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.