




Recombinant Human Serotransferrin (TF) (HOLO) (Active)

Product Code	CSB-MP023412HU-WD1
Abbreviation	Recombinant Human TF protein (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P02787
Form	Lyophilized
Storage Buffer	Lyophilized from a 0.2 µm filtered solution in sterile water.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Biological Activity	Measured in a serum-free cell proliferation assay using MDCK canine kidney epithelial cells. The ED50 for this effect is 0.5-2 µg/mL.
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	VPDKTVRWCAVSEHEATKQCQSFDRDHMKSVIPSDGPSVACVKKASYLDCIRAI ANEADAVTL DAGLVYDAYLAPNNLKPVVAEFYGSKEDPQTFYYAVAVVKKDSG FQMNQLRGKKSCHTGLGRSAGWNIPIGLLYCDLPEPRKPLEKAVANFFSGSC APCADGTD FFPQLCQLCPGCGCSTLNQYFGYSGAFKCLKDGAGDVAFVKHSTI FENLANKARDQYELLCLDNTRKPVDEYKDCHLAQVPSHTVVARSMGGKEDL IWELLNQAQEHFGKDKSKEFQLFSSPHGKDLLFKDSA HGF LKVP PRMDAKMY LGYEYVTAIRNLREGTCPEAPTDECKPVKWCALSHHERLKCDEWSVNSVGKIE CVSAETTEDCIAKIMNGEADAMSLDGGFVYIAGKCGLVPVLAENYNKSDNCED TPEAGYFAIAVVKKSASDLTWDNLKGKKSCHTAVGRTAGWNIPMGLLYNKINH CRFDEFFSEGCAPGSKKDSSLCKLCMGSGNLNCEPNNKEGYGYTGAFRCLV EKGDVAFVKHQTVPQNTGGKNPDPWAKNLNEKDYELLCLDGTRKPV E EYAN CHLARAPNHAVVTRKDKEACVHKILRQQQHLFGSNVTDCSGNFCLFRSETKD LLFRDDTVCLAKLHDRNTYEKYL GEEYVKA VGNLRK CSTSSLLEACTFRRP
Source	Mammalian cell
Target Names	Transferrin
Expression Region	20-698aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag free
Mol. Weight	75.2 kDa
Protein Length	Full Length of Mature Protein



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Endotoxin

≤0.5 EU/mg by the LAL method

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.