




Recombinant Mycobacterium tuberculosis NAD(P)H dehydrogenase (quinone) (IpdA)

Product Code	CSB-EP6072FSG
Abbreviation	Recombinant Mycobacterium tuberculosis IpdA protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P9WHH7
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Mycobacterium tuberculosis NAD(P)H dehydrogenase (quinone) (IpdA)
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MVTRIVILGGGPAGYEAALVAATSHPETTQVTVIDCDGIGGA AVLDDCVPSKTFI ASTGLRTELRRAPHLGFHIDFDDAKISLPQIHARVKTLAAAQSADITAQLLSMGV QVIAGRGELIDSTPGLARHRIKATAADGSTSEHEADVVLVATGASPRILPSAQP DGERILTWRQLYDLDPDHLIVVGSGVTGAEFVDAYTELGVPTTVASQDHV LPYEDADAALVLEESFAERGVR LFKNARAASVTRTGAGVLVTMTDGRTVEGS HALMTIGSV PNTSGLGLERVGIQLGRGNLYTVDRVSRTLATGIYAAGDCTGLLP LASVAAMQGRIAMYHALGEGVSP IRLRTVAATVTRPEIAAVGV PQSVIDAGSV AARTIMLPLRTNARAKMSEMRHGFVKIFCRRSTGVVIGGVVVAPIASELILPIAV AVQNRITVNELAQTLAVYPSLSGSITEAARRLMAHDDL DCTAAQDAAEQLALVP HHLPTSN
Research Area	Others
Source	E.coli
Target Names	IpdA
Expression Region	1-493aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	58.4 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.