




Whitlow linker Monoclonal Antibody, FITC conjugated

Product Code	CSB-MA234468A0m-C
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Whitlow linker-KLH
Species Reactivity	Homo sapiens (Human)
Tested Applications	FC; Recommended dilution: FC:1:50-1:200
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG2b
Clonality	Monoclonal
Product Type	Monoclonal Antibody
Research Area	Others
Target Names	Whitlow linker
Clone No.	18A3G5

Image



Untransfected HEK-293F cells (green line) and transfected scFv-based Anti-CD19 CAR containing a Whitlow/218 linker stable cells (red line) were stained with anti-Whitlow linker FITC-conjugated antibody (2µg/1*10⁶cells).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.