



G4S linker Monoclonal Antibody, FITC conjugated

Product Code	CSB-MA186062I3m-C
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	G4S-BSA
Species Reactivity	Homo sapiens (Human)
Tested Applications	FC; Recommended dilution: FC:1:50-1:200
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG2a
Clonality	Monoclonal
Product Type	Monoclonal Antibody
Research Area	Others
Target Names	G4S linker
Clone No.	14A10E7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.