



HEK293/Human CLDN9-eGFP Stable Cell Line

Product Code	CSB-SC005511HU
Target Name	Adherent
Storage	Frozen in liquid nitrogen or -80°C
Uniprot No.	O95484
Tested Applications	Binding assay by FACS and cell based ELISA.
Storage Buffer	35% FBS + 55% DMEM + 10% DMSO
Resistance	Puromycin
Intended Use	DMEM + 10% FBS + 0.75 µg/mL Puromycin
Tag Info	C-terminal eGFP-tagged

Image



Untransfected HEK-293 cells (green line) and transfected Human CLDN9 HEK-293 stable cells (red line) were stained with anti-CLDN6/9 antibody(CSB-RA005508MA1HU) (2µg/1*10⁶cells), washed and then followed by APC-conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.



Untransfected HEK-293 cells (green line) and transfected Human CLDN9 HEK-293 stable cells (red line) were stained with anti-CLDN6/9 antibody(CSB-RA005508MA3HU) (2µg/1*10⁶cells), washed and then followed by APC-conjugated anti-Mouse IgG Fc antibody and analyzed with flow cytometry.