




HEK293T/Human KIR3DL2 Stable Cell Line

Product Code	CSB-SC012365HU1
Target Name	Adherent
Storage	Frozen in liquid nitrogen or -80°C
Uniprot No.	P43630
Tested Applications	Binding assay by FACS and cell based ELISA.
Storage Buffer	35% FBS + 55% DMEM + 10% DMSO
Resistance	Puromycin
Intended Use	DMEM + 10% FBS + 0.5 µg/ml Puromycin
Tag Info	N-terminal Flag-tagged

Image



Untransfected HEK-293T cells (green line) and transfected Human KIR3DL2 HEK-293T stable cells (red line) were stained with anti-KIR3DL2 antibody (CSB-RA012365MA1HU) (2 µg/1*10⁶ cells), washed and then followed by FITC conjugated anti-Human IgG Fc antibody and analyzed with flow cytometry.