



Palmitoyl Pentapeptide-4

Product Code	CSB-DT1698
Abbreviation	214047-00-4
Storage	Store at -20°C upon receipt, aliquoting is necessary for multiple use. Avoid repeated freeze-thaw cycles. Short term: 4 °C, one week from the date of receipt. Long term: -20 °C to -80 °C, 6 months from the date of receipt.
Form	Lyophilized powder
Purity	98%
Sequence	Pal-Lys-Thr-Thr-Lys-Ser
Lead Time	1-2 weeks
Resistance	Pal-KTTKS
Notes	The shelf life is related to many factors, storage state, storage temperature and the stability of the product itself. The above shelf lives are concluded by our current storage experience. The specific shelf life will depend on specific product. Please use as soon as possible after receiving. Repeated freezing and thawing is not recommended.
Mol. Weight	802.065 Da
Accession NO.	C39H75N7O10