



Palmitoyl Tripeptide-8

Product Code	CSB-DT1701
Abbreviation	936544-53-5
Storage	Store at -20°C upon receipt, aliquoting is necessary for multiple use. Avoid repeated freeze-thaw cycles. Short term: 4 °C, one week from the date of receipt. Long term: -20 °C to -80 °C, 6 months from the date of receipt.
Form	Lyophilized powder
Purity	98%
Sequence	Pal-His-D-Phe-Arg-NH ₂
Lead Time	1-2 weeks
Resistance	Pal-HDFR-NH ₂
Notes	The shelf life is related to many factors, storage state, storage temperature and the stability of the product itself. The above shelf lives are concluded by our current storage experience. The specific shelf life will depend on specific product. Please use as soon as possible after receiving. Repeated freezing and thawing is not recommended.
Mol. Weight	695.949 Da
Accession NO.	C37H61N9O4