



a-MSH

Product Code	CSB-DT1738
Abbreviation	581-05-5
Storage	Store at -20°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles. Short term: 4 °C, one week from the date of receipt. Long term: -20 °C to -80 °C, 6 months from the date of receipt.
Form	Lyophilized powder
Purity	98%
Sequence	Ac-Ser-Tyr-Ser-Met-Glu-His-Phe-Arg-Trp-Gly-Lys-Pro-Val-NH2
Lead Time	1-2 weeks
Resistance	Ac-SYSMEHFRWGKPV-NH2
Notes	The shelf life is related to many factors, storage state, storage temperature and the stability of the product itself. The above shelf lives are concluded by our current storage experience. The specific shelf life will depend on specific product. Please use as soon as possible after receiving. Repeated freezing and thawing is not recommended.
Mol. Weight	1664.9 Da
Accession NO.	C77H109N21O19S1