



N-(3-(3-hydroxyprop-1-yn-1-yl)-5-(trifluoromethyl)pyridin-2-yl)pivalamide

Product Code	CSB-DT4935
Abbreviation	1346447-01-5
Storage	Store at -20°C upon receipt, aliquoting is necessary for mutiple use. Avoid repeated freeze-thaw cycles. Short term: 4 °C, one week from the date of receipt. Long term: -20 °C to -80 °C, 6 months from the date of receipt.
Form	Lyophilized powder
Purity	95%
Lead Time	1-2 weeks
Notes	The shelf life is related to many factors, storage state, storage temperature and the stability of the peptide itself. The above shelf lives are concluded by our current storage experience. The specific shelf life will depend on specific product. Please use as soon as possible after receiving. Repeated freezing and thawing is not recommended.
Accession NO.	C14H15F3N2O2